

STEP BY STEP INSTRUCTONS
I.V. MINI BAG WITH HOHN/SUBCLAVIAN

1. GATHER SUPPIES: MEDICATIONS, SODIUM CHLORIDE AND HEPARIN FLUSHES, AND ANTIMICROBIAL SOAP.
2. WASH HANDS WITH ANTIMICROBIAL SOAP.
3. OBTAIN NEEDED SODIUM CHLORIDE AND HEPARIN FLUSHES AND REMOVE FROM WRAPPERS.
4. OBTAIN IV MINI-BAG:_____ SQUEEZE DRIP CHAMBER UNTIL APPROX. ½ FULL. OPEN THE ROLLER CLAMP AND LET THE TUBING BECOME FILLED WITH FLUID. ONCE THE FLUID HAS REACHED THE END OF TUBING, CLOSE THE ROLLER CLAMP.
5. REMOVE BLUE CAP FROM END OF CATHETER AND UNCLAMP CATHETER.
6. OBTAIN SODIUM CHLORIDE, REMOVE END CAP AND FLUSH OUT EXCESS AIR.
7. TWIST SYRINGE ONTO END OF CATHETER, FLUSH AND DISCARD INTO HOME TRASH RECEPTACLE.
8. REMOVE CAP FROM THE END OF THE ANTIBIOTIC TUBING AND TWIST ONTO THE END OF THE CATHETER.
9. OPEN THE ROLLER CLAMP ON THE TUBING ENOUGH TO ALLOW THE MEDICATION TO INFUSE AT ABOUT 1 DROP PER SEC. MEDICATION SHOULD INFUSE IN ABOUT _____ MIN/ HRS.
10. ONCE THE MEDICATION IS INFUSED, CLOSE THE ROLLER CLAMP. REMOVE ANTIBIOTIC TUBING FROM THE END OF THE CATHETER AND DISCARD.
11. OBTAIN SODIUM CHLORIDE, REMOVE END CAP AND FLUSH OUT EXCESS AIR. TWIST SYRINGE ONTO END OF CATHETER, FLUSH AND DISCARD.
12. OBTAIN HEPARIN, REMOVE END CAP AND FLUSH OUT EXCESS AIR. TWIST SYRINGE ONTO END OF CATHETER, FLUSH AND DISCARD.
13. PLACE BLUE CAP ONTO END OF CATHETER AND CLOSE CLAMP.
14. DISCARD ALL USED SUPPLIES INTO HOME TRASH RECEPTACLE.
15. FLUSH UNUSED CATHETER LUMEN(S) WITH HEPARIN ONCE DAILY.